OXYBOX

SYRGASLÅDA A1

MANUAL

INNEHÅLL

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SAFETY REGULATIONS

**Oxygen, together with oil and grease provides increased fire risk and can cause an explosion!
All details **MUST** be FAT- and OIL-free!**

This equipment delivers 100% oxygen and should be used with caution in the present of a fire!

Demand valves purge function (demand button), if available, may only be used by trained personnel and only on adult patients.

Repairs, maintenance, modification or other alternations may only be performed by technical personnel with special training.

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Pos	Art.nr	Title
1	60-100	Oxygenregulator S.O.S
2	60-10102	Pocketmask
3	60-1011	Oxygenmask/silicon
4	60-1013	Constant flow mask
5	60-1045	LP-hose O ₂ / 2,5m
6	60-6210	5L steel cylinder w.o / valve
7	60-6017	Velcro bracket
8	60-6020	Case interior 821
9	60-60401	Case gland
10	60-6215	Valve M90
11	60-6345	Demandvalve S.O.S
12	68-003092	Box 309 Yellow
13	68-008212	Ultracase 821 Yellow

TECHNICAL DATA

Case

Weight:	13 kg
Case material:	ABS-plastic
Temp range:	-30 to +60 °C
Dimensions:	W=47 cm, D=37 cm, H=23 cm

Reducing valve

Cylinder pressure:	5-220 bar
O ₂ quick connection outlet:	approx 5 bar (max 85 l/min)
Low pressure flow outlet:	0-25 l/min in 12 steps
Pressure gauge:	Scale: 0-315 bar

Oxygen cylinder 5 l

Filling pressure:	200 bar
Weight-filled:	9,5 kg
Oxygen mass:	1000 l
Treatment time:	65 min (at 15 l/min)

Demandvalve

Supply pressure:	2,8-6,2 bar
Flow range:	0-160 l/min
Temp range :	-30 to + 50 °C
Dead space:	8 ml (ex. mask)

FUNCTIONAL CHECK


ALWAYS check that the reduction valves low pressure selector is in position "0" (All other settings means that oxygen is lost through the low pressure outlet)

Check the equipment by:

- a) Connect the white oxygen hose to the quick connection on the reduction valve, linking the reduction valve together with the demand valve.
- b) Carefully open the valve on the cylinder.
Pressure gauge will then show a pressure between 5 and 200 bar depending on the degree of filling.
- c) Check audio-visual that there are no leaks.
No leaks may exist.
- d) Test the oxygen therapy equipment by connect the oxygen mask to the demand valve, hold the mask against your own face and take 3 deep breaths. If the system works without any faults, clean the oxygen mask and put the oxygen therapy equipment in "standby mode" or unpack the equipment for transport.
- e) After a complet function check, put the oxygen therapy equipment in "standby mode" by closing the valve on the cylinder. Residual pressure in the system can be emptied through the low pressure outlet on the reduction valve. Open the switch and let all remaing oxygen out through the low pressure outlet. Then close the flow switch by setting it in "0" position.

QUICK GUIDE

(from "standby mode" to treatment of patient)

OXYGEN TREATMENT OF PATIENT

- 1) Prepare oxygen therapy equipment according to previously specified instructions.
- 2) Check the available amount of oxygen.
- 3) Make sure the patient's airways is free.
- 4) Explain to the patient
 - He / she shall receive oxygen through a mask
 - The oxygen is part of the treatment
 - To breathe normally through the mask.
 - He / she should notify if something feels strange.

Place the breathing mask tightly against the patient's face, over the nose and mouth. Let the patient breathe normally at his/her own pace if this is deemed sufficient. If not, the patient requires breathing assistance / artificial respiration.

See part 8 for more information for patients with inadequate breathing or respiratory arrest.

- 5) The patient must not be left unattended and remember to now and then, check the function, for oxygen leaks, and that no obstacles exist in the airways in the form of secretions and vomiting.
- 6) Ensure that good ventilation occurs in the area around the patient so that no leakage of oxygen and the patient's exhaled air accumulate in loose clothing, blankets, etc., which increase the risk of fire!

- 7) **No fire or sparking objects may be present in the vicinity of the patient!**
- 8) **Artificial oxygen therapy can only be performed by trained personnel with patients having reduced and / or inadequate breathing by using the purge button on the demand valve, if one is available.**

One patient with no or inadequate breathing must always first be suspected to have obstructed airway. It is therefore very important to immediately establish a free airway.

For patients with respiratory arrest, it is vitally important that artificial respiration starts immediately. Please start mouth-to-mouth breathing while waiting for the Oxyboxen to be prepared.

Pocket mask can be used for mouth-to-mouth treatment, it has a built-in filter and protects the user against health hazards.

Note that misuse of the purge button while treating a patient highly increase the risk of damaging the respiratory tract due to overpressure!